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SILVERLON**PERFORMANCE. POWER. SPEED.****Device Description**

Silverlon® Acute Burn Glove Dressings are a sterile, single layer, non-adherent, antimicrobial barrier dressing in the form of a glove with open finger, which is designed to fit comfortably over the hand, covering the entire skin surface of the hand.

Available as:

Code	Size
ABG-01XS	Adult Extra Small
ABG-01S	Adult Small
ABG-01M	Adult Medium
ABG-01L	Adult Large
ABG-01XL	Adult Extra Large
ABG-01XXL	Adult Extra Extra Large

Silverlon® Acute Burn Glove Dressing contains silver to provide an ancillary, local antimicrobial effect to reduce the risk of wound infection by organisms sensitive to silver. Ionic silver is a broad spectrum antimicrobial, which has been shown to be effective against a range of pathogens associated with wound infection. Silverlon® dressings have been tested in vitro and found effective against microorganisms such as: Enterococcus faecalis, Staphylococcus aureus, Pseudomonas aeruginosa, Acinetobacter Baumannii, Enterobacter cloacae, Staphylococcus epidermidis and Klebsiella pneumoniae.

In vitro testing has shown Silverlon® dressings to be MRI compatible.

Indications

Silverlon® Acute Burn Glove Dressing is a professional use wound dressing intended for local management of the following types of hand wounds:

- 1st and 2nd degree burns
- Donor and graft site wounds
- Traumatic and surgical wounds
- Infected wounds*

*Silverlon® Acute Burn Glove Dressing may be used on those infected wound types listed above. Where used on infected wounds, the infection should be treated as per local clinical protocol.

Silverlon® Acute Burn Glove Dressing is indicated for a total duration of contact (i.e. where treatment involves consecutive application of individual dressings) of up to 30 days.

Warnings

- Do not use past expiration date on the product packaging.
- Do not use if pouch is damaged or open.
- Clinicians / Healthcare Professionals should be aware that there are very limited data on prolonged and repeated use of silver-containing dressings, particularly in children and neonates.
- Do not use petroleum-based ointments or creams under Silverlon® Acute Burn Glove Dressing.
- Do not moisten Silverlon® Acute Burn Glove Dressing with hydrogen peroxide or povidone iodine.
- Silverlon® Acute Burn Glove Dressings are intended for single use only and should not be re-used. Re-use may adversely affect the dressing’s performance characteristics and also presents a risk of infection to the patient.
- Do not remove Silverlon® Acute Burn Glove Dressing from the hand while dry or if the dressing sticks when trying to remove it. Follow the instructions below for dressing removal.
- Total duration of use of Silverlon® Acute Burn Glove Dressing (i.e. where treatment involves consecutive application of individual dressings) should not exceed 30 days.

LC-CE-IFU-ABG-02; Rev B**Emergo Europe, Westervoortsedijk 60; 6827 AT Arnhem The Netherlands Tel. +(31) (0) 70 345-8570****Argentum Medical, LLC 2571 Kaneville Ct Geneva, IL 60134 USA Tel. (+1) 888.551.0188 Fax. (+1) 888.558.9923 www.silverlon.com****SILVERLON****PERFORMANCE. POWER. SPEED.****Device Description**

Silverlon® Acute Burn Glove Dressings are a sterile, single layer, non-adherent, antimicrobial barrier dressing in the form of a glove with open finger, which is designed to fit comfortably over the hand, covering the entire skin surface of the hand.

Available as:

Code	Size
ABG-01XS	Adult Extra Small
ABG-01S	Adult Small
ABG-01M	Adult Medium
ABG-01L	Adult Large
ABG-01XL	Adult Extra Large
ABG-01XXL	Adult Extra Extra Large

Silverlon® Acute Burn Glove Dressing contains silver to provide an ancillary, local antimicrobial effect to reduce the risk of wound infection by organisms sensitive to silver. Ionic silver is a broad spectrum antimicrobial, which has been shown to be effective against a range of pathogens associated with wound infection. Silverlon® dressings have been tested in vitro and found effective against microorganisms such as: Enterococcus faecalis, Staphylococcus aureus, Pseudomonas aeruginosa, Acinetobacter Baumannii, Enterobacter cloacae, Staphylococcus epidermidis and Klebsiella pneumoniae.

In vitro testing has shown Silverlon® dressings to be MRI compatible.

Indications

Silverlon® Acute Burn Glove Dressing is a professional use wound dressing intended for local management of the following types of hand wounds:

- 1st and 2nd degree burns
- Donor and graft site wounds
- Traumatic and surgical wounds
- Infected wounds*

*Silverlon® Acute Burn Glove Dressing may be used on those infected wound types listed above. Where used on infected wounds, the infection should be treated as per local clinical protocol.

Silverlon® Acute Burn Glove Dressing is indicated for a total duration of contact (i.e. where treatment involves consecutive application of individual dressings) of up to 30 days.

Warnings


- Do not use past expiration date on the product packaging.
- Do not use if pouch is damaged or open.
- Clinicians / Healthcare Professionals should be aware that there are very limited data on prolonged and repeated use of silver-containing dressings, particularly in children and neonates.
- Do not use petroleum-based ointments or creams under Silverlon® Acute Burn Glove Dressing.
- Do not moisten Silverlon® Acute Burn Glove Dressing with hydrogen peroxide or povidone iodine.
- Silverlon® Acute Burn Glove Dressings are intended for single use only and should not be re-used. Re-use may adversely affect the dressing’s performance characteristics and also presents a risk of infection to the patient.
- Do not remove Silverlon® Acute Burn Glove Dressing from the hand while dry or if the dressing sticks when trying to remove it. Follow the instructions below for dressing removal.
- Total duration of use of Silverlon® Acute Burn Glove Dressing (i.e. where treatment involves consecutive application of individual dressings) should not exceed 30 days.

Contraindications

- Avoid using Silverlon® Acute Burn Glove Dressing on patients with known sensitivity to silver or nylon.

SILVERLON**PERFORMANCE. POWER. SPEED.**

Silverlon® Acute Burn Glove Dressings

Publication Date: 5/2023**Emergo Europe, Westervoortsedijk 60; 6827 AT Arnhem The Netherlands Tel. +(31) (0) 70 345-8570****Argentum Medical, LLC 2571 Kaneville Ct Geneva, IL 60134 USA Tel. (+1) 888.551.0188 Fax. (+1) 888.558.9923 www.silverlon.com****LC-CE-IFU-ABG-02; Rev B****Silverlon® Acute Burn Glove Dressings****English**

- Silverlon® Acute Burn Glove Dressing is not intended for use on 3rd degree burns.

Adverse Reactions

Silverlon® Acute Burn Glove Dressing has been subjected to independent standard in vitro and in vivo biocompatibility tests, including cytotoxicity, sensitization, and intracutaneous reactivity. All tests were performed in accordance with International Standard Organization (ISO) 10993 Standard Series for Biological Evaluation of Medical Devices. The results of the studies indicated that Silverlon® Acute Burn Glove Dressing is safe for its intended use.

Instructions for Use

- Cleanse the wound with sterile water, distilled water, or normal saline removing necrotic debris or eschar as needed per local protocol.
- Select the Silverlon® Acute Burn Glove Dressing size that best covers the hand and makes light contact with the skin.
- The glove fingers may be trimmed to fit (or folded over and the end secured) per local protocols.
- Activate Silverlon® Acute Burn Glove Dressing by thoroughly moistening it with sterile water, distilled water or normal saline.
- Gently slide the hand into the glove and secure the glove in position with an outer secondary dressing (i.e. gauze wrap, tubular stretch knit bandage or other compressive wrap dressing) per local protocol.
 - It is recommended to use a moisture-donating outer dressing such as hydrocolloid or pre-moistened foam or gauze.
 - Wrap and secure the outer secondary dressing as per local protocol.
- Silverlon® Acute Burn Glove Dressing may be worn while bathing / showering, although the outer secondary dressing should be removed first.
- Silverlon® Acute Burn Glove Dressing is intended for use for up to 7 days.
- To remove Silverlon® Acute Burn Glove Dressing to enable wound inspection and clinical attention, first remove and discard the outer secondary dressing, then soak the hand encased in Silverlon® Acute Burn Dressing in room temperature sterile water, distilled water, or normal saline for 5 -15 minutes.
 - Next, gently remove Silverlon® Acute Burn Dressing.
 - If the dressing sticks to the skin or wound during removal, repeat soaking until it can be easily removed.

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