SILVERLON[®]
safe, strong, simple™

Apósitos tipo guante para quemaduras graves Silverlon®

Spanish

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English

Device Description

Silverlon® Acute Burn Glove Dressings are a sterile, single layer, non-adherent, antimicrobial barrier dressing in the form of a glove with open finger, which is designed to fit comfortably over the hand, covering the entire skin surface of the hand.

Available as:

Code	Size
ABG-01XS	Adult Extra Small
ABG-01S	Adult Small
ABG-01M	Adult Medium
ABG-01L	Adult Large
ABG-01XL	Adult Extra Large
ABG-01XXL	Adult Extra Extra Large

Silverlon® Acute Burn Glove Dressing contains silver to provide an ancillary, local antimicrobial effect to reduce the risk of wound infection by organisms sensitive to silver. Ionic silver is a broad spectrum antimicrobial, which has been shown to be effective against a range of pathogens associated with wound infection. Silverlon® dressings have been tested in vitro and found effective against micro organisms such as: Enterococcus faecalis, Staphylococcus aureus, Pseudomonas aeruginosa, Acinetobacter Baumannii, Enterobacter cloacae, Staphylococcus epidermidis and Klebsiella pneumoniae.

In vitro testing has shown Silverlon® dressings to be MRI compatible.

Indications

Silverlon® Acute Burn Glove Dressing is a professional use wound dressing intended for local management of the following types of hand wounds:

- 1st and 2nd degree burns
- Donor and graft site wounds
- Traumatic and surgical wounds
- Infected wounds*

*Silverlon® Acute Burn Glove Dressing may be used on those infected wound types listed above. Where used on infected wounds, the infection should be treated as per local clinical protocol.

Silverlon® Acute Burn Glove Dressing is indicated for a total duration of contact (i.e. where treatment involves consecutive application of individual dressings) of up to 30 days.


Warnings

- Do not use past expiration date on the product packaging.
- Do not use if pouch is damaged or open.
- Clinicians / Healthcare Professionals should be aware that there are very limited data on prolonged and repeated use of silver-containing dressings, particularly in children and neonates.
- Do not use petroleum-based ointments or creams under Silverlon® Acute Burn Glove Dressing.
- Do not moisten Silverlon® Acute Burn Glove Dressing with hydrogen peroxide or povidone iodine.
- Silverlon® Acute Burn Glove Dressings are intended for single use only and should not be re-used. Re-use may adversely affect the dressing's performance characteristics and also presents a risk of infection to the patient.
- Do not remove Silverlon® Acute Burn Glove Dressing from the hand while dry or if the dressing sticks when trying to remove it. Follow the instructions below for dressing removal.
- Total duration of use of Silverlon® Acute Burn Glove Dressing (i.e. where treatment involves consecutive application of individual dressings) should not exceed 30 days.

SILVERLON[®]

safe, strong, simple™

Silverlon® Acute Burn Glove Dressings

Emergo Europe
Molenstraat 15
2513 BH The Hague, The Netherlands
Tel: +(31) (0) 70 345-8570
Fax: +(31) (0) 70 346-7299

Argentum Medical
2571 Kaneville Ct
Geneva, IL 60134 USA
Tel. (+1) 888.551.0188
Fax. (+1) 888.558.9923
www.silverlon.com
LC-CE-IFU-ABG-Rev 03