**SILVERLON®**  
safe, strong, simple

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**SILVERLON®**  
safe, strong, simple

## Silverlon® Lifesaver® IV/Catheter Dressings

**English**

### Device Description

Silverlon® Lifesaver® IV/ Catheter Dressings are sterile, multi-layer, non-adherent, absorbent, antimicrobial barrier dressings, intended for topical use in wound management, to reduce risk of wound infection.

Available as:

<b>Code</b>	<b>Pad Size</b>
Lifesaver™ Ag IV/CD-15	2,5 cm circle with 1.5 mm center hole
Lifesaver™ Ag IV/CD-40	2,5 cm circle with 4.0 mm center hole
Lifesaver™ Ag IV/CD-70	2,5 cm circle with 7.0 mm center hole

Silverlon® Lifesaver® IV/ Catheter Dressings contain silver to provide an ancillary, local antimicrobial effect to reduce the risk of wound infection by organisms sensitive to silver. Ionic silver is a broad spectrum antimicrobial, which has been shown to be effective against a range of pathogens associated with wound infection. Silverlon® dressings have been tested in vitro and found effective against micro organisms such as: Enterococcus faecalis, Staphylococcus aureus, Pseudomonas aeruginosa, Acinetobacter Baumannii, Enterobacter cloacae, Staphylococcus epidermidis, and Klebsiella pneumoniae.

In vitro testing has shown Silverlon® dressings to be MRI compatible.

### Indications

Silverlon® Lifesaver® IV/ Catheter Dressings are professional use wound dressings indicated for local management of:

- vascular access sites (including around IV catheters, central venous lines, arterial catheters, dialysis catheters);
- drainage tube access sites;
- external orthopedic pin sites, and
- epidural catheter access sites.

Silverlon® Lifesaver® IV/ Catheter Dressings are indicated for a total duration of contact (i.e. where treatment involves consecutive application of individual dressings) of up to 30 days.

Silverlon® Lifesaver® IV/Catheter Dressings may be used on infected wounds, under medical supervision, in conjunction with systemic antimicrobial therapy.

### Warnings

- Do not use past expiration date on the product packaging.
- Do not use if pouch is damaged or open.
- Clinicians / Healthcare Professionals should be aware that there are very limited data on prolonged and repeated use of silver-containing dressings, particularly in children and neonates.
- Silverlon® dressings should not be used during pregnancy or lactation because of the lack of data concerning the effects of silver on reproduction.
- Do not use petroleum-based ointments or creams under Silverlon® Lifesaver® IV/ Catheter Dressing.
- Do not soak Silverlon® Lifesaver® IV/Catheter Dressing in hydrogen peroxide, or povidone iodine.
- Silverlon® Lifesaver® IV/Catheter Dressings are intended for single use only and should not be re-used. Re-use may adversely affect the dressing’s performance characteristics and also presents a risk of infection to the patient.
- Total duration of use of Silverlon® Lifesaver® IV/Catheter Dressing (i.e. where treatment involves consecutive application of individual dressings) should not exceed 30 days.

### Contraindications

- Avoid using Silverlon® Lifesaver® IV/Catheter Dressings on patients with known sensitivity to silver or nylon.

**SILVERLON®**

safe, strong, simple™

# Silverlon® Lifesaver® IV/Catheter Dressings

Emergo Europe  
Molenstraat 15  
2513 BH The Hague, The Netherlands  
Tel: +(31) (0) 70 345-8570  
Fax: +(31) (0) 70 346-7299

Argentum Medical  
2571 Kaneville Ct  
Geneva, IL 60134 USA  
Tel. (+1) 888.551.0188  
Fax. (+1) 888.558.9923  
www.silverlon.com  
LC-CE-IFU-IVCD Rev. 02

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