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LC-CE-IFU-NPS-FLX Rev 02

SILVERLON
PERFORMANCE. POWER. SPEED.

Silverlon® flex Wound Contact Negative Pressure Dressings

English

Device Description

Silverlon® flex Wound Contact Negative Pressure Dressings are sterile, single layer, non-adherent, antimicrobial barrier dressings, intended for topical use in wound management and in negative pressure wound therapy (NPWT), to reduce risk of wound infection.

Available as:

Silverlon® flex Wound Contact Negative Pressure Dressings

Code	Size(cm)	Code	Size (cm)
WCD22-NPS	5 x 5 cm	WCD48-NPS	10 x 20 cm
WCD44-NPS	10 x 10 cm	WCD88-NPS	20 x 20 cm

Silverlon® flex Wound Contact Negative Pressure Dressings contain silver to provide an ancillary, local antimicrobial effect to reduce the risk of wound infection by organisms sensitive to silver. Ionic silver is a broad spectrum antimicrobial, which has been shown to be effective against a range of pathogens associated with wound infection.

Silverlon® dressings have been tested in vitro and found effective against microorganisms such as: Enterococcus faecalis, Staphylococcus aureus, Pseudomonas aeruginosa, Acinetobacter Baumannii, Enterobacter cloacae, Staphylococcus epidermidis and Klebsiella pneumoniae. Silverlon® flex Wound Contact Negative Pressure Dressings may be used with all types of negative pressure wound therapy porous foam/gauze dressings and equipment.

In vitro testing has shown Silverlon® flex dressings to be MRI compatible.

Indications

Silverlon® flex Wound Contact Negative Pressure Dressings are professional use wound dressings indicated for wound management and may also be used in conjunction with negative pressure wound therapy (NPWT), for local management of:

- 1st and 2nd degree burns
- Donor and graft site wounds
- Diabetic, pressure and venous ulcers
- Traumatic and surgical wounds
- Infected wounds*

* Silverlon® flex Wound Contact Negative Pressure Dressings may be used on those infected wounds listed above. Where used on infected wounds, the infection should be treated as per local clinical protocol.

Silverlon® flex Wound Contact Negative Pressure Dressings are indicated for a total duration of contact (i.e. where treatment involves consecutive application of individual dressings) of up to 30 days.

Warnings

- Do not use past expiration date on the product packaging.
- Do not use if pouch is damaged or open.
- Clinicians / Healthcare Professionals should be aware that there are very limited data on prolonged and repeated use of silver-containing dressings, particularly in children and neonates.
- Silverlon® flex Wound Contact Negative Pressure Dressings should not be used during pregnancy or lactation because of the lack of data concerning the effects of silver on reproduction.
- Silverlon® flex Wound Contact Negative Pressure Dressings are intended for single use only and should not be re-used. Re-use may adversely affect the dressing's performance characteristics and also presents a risk of infection to the patient.
- Do not use petroleum-based ointments or creams under Silverlon® flex Wound Contact Negative Pressure Dressings.
- Do not moisten Silverlon® flex Wound Contact Negative Pressure Dressings with hydrogen peroxide or povidone iodine.
- Total duration of use of Silverlon® flex Wound Contact Negative Pressure Dressings (i.e. where treatment involves consecutive application of individual dressings) should not exceed 30 days.
- Silverlon® flex Wound Contact Negative Pressure Dressings are not intended for direct application to exposed underlying structures (i.e. blood vessels, bone, tendons, nerves).

Compressas Silverlon® flex Wound Contact Negative Pressure

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