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
ENGLISH**SILVERLON****PERFORMANCE. POWER. SPEED.**

Silverlon® Flex Wound Contact Negative Pressure Dressings

**EC****REP**

Emergo Europe, Westervoortsedijk 60; 6827 AT Arnhem The Netherlands Tel. +31) (0) 70 345-8570

LC-CE-IFU-NPS-FLX-03; Rev B**SILVERLON****PERFORMANCE. POWER. SPEED.****Device Description**

Silverlon® flex Wound Contact Negative Pressure Dressings are sterile, single layer, non-adherent, antimicrobial barrier dressings, intended for topical use in wound management and in negative pressure wound therapy (NPWT), to reduce risk of wound infection.

Available as:

Silverlon® flex Wound Contact Negative Pressure Dressings

Code	Size	Code	Size
WCD22-NPS(-FLX)	5 x 5 cm	WCD48-NPS(-FLX)	10 x 20 cm
WCD44-NPS(-FLX)	10 x 10 cm	WCD88-NPS(-FLX)	20 x 20 cm

Silverlon® flex Wound Contact Negative Pressure Dressings contain silver to provide an ancillary, local antimicrobial effect to reduce the risk of wound infection by organisms sensitive to silver. Ionic silver is a broad spectrum antimicrobial, which has been shown to be effective against a range of pathogens associated with wound infection. Silverlon® dressings have been tested in vitro and found effective against microorganisms such as: Enterococcus faecalis, Staphylococcus aureus, Pseudomonas aeruginosa, Acinetobacter Baumannii, Enterobacter cloacae, Staphylococcus epidermidis and Klebsiella pneumoniae. Silverlon® flex Wound Contact Negative Pressure Dressings may be used with all types of negative pressure wound therapy porous foam/gauze dressings and equipment. In vitro testing has shown Silverlon® flex dressings to be MRI compatible.

Indications

Silverlon® flex Wound Contact Negative Pressure Dressings are professional use wound dressings indicated for wound management and may also be used in conjunction with negative pressure wound therapy (NPWT), for local management of:

- 1st and 2nd degree burns
- Donor and graft site wounds
- Diabetic, pressure and venous ulcers
- Traumatic and surgical wounds
- Infected wounds*

* Silverlon® flex Wound Contact Negative Pressure Dressings may be used on those infected wounds listed above. Where used on infected wounds, the infection should be treated as per local clinical protocol.

Silverlon® flex Wound Contact Negative Pressure Dressings are indicated for a total duration of contact (i.e. where treatment involves consecutive application of individual dressings) of up to 30 days.

Warnings

- Do not use past expiration date on the product packaging.
- Do not use if pouch is damaged or open.
- Clinicians / Healthcare Professionals should be aware that there are very limited data on prolonged and repeated use of silver-containing dressings, particularly in children and neonates.
- Silverlon® flex Wound Contact Negative Pressure Dressings should not be used during pregnancy or lactation because of the lack of data concerning the effects of silver on reproduction.
- Silverlon® flex Wound Contact Negative Pressure Dressings are intended for single use only and should not be re-used. Re-use may adversely affect the dressing’s performance characteristics and also presents a risk of infection to the patient.
- Do not use petroleum-based ointments or creams under Silverlon® flex Wound Contact Negative Pressure Dressings.
- Do not moisten Silverlon® flex Wound Contact Negative Pressure Dressings with hydrogen peroxide or povidone iodine.
- Total duration of use of Silverlon® flex Wound Contact

© 2023 Emergo Europe. All rights reserved. Silverlon® and Silverlon® flex are registered trademarks of Emergo Europe. All other trademarks are the property of their respective owners.

English

Negative Pressure Dressings (i.e. where treatment involves consecutive application of individual dressings) should not exceed 30 days.

- Silverlon® flex Wound Contact Negative Pressure Dressings are not intended for direct application to exposed underlying structures (i.e. blood vessels, bone, tendons, nerves).
- When using Silverlon® flex Wound Contact Negative Pressure Dressings in conjunction with commercially available negative pressure wound therapy systems, users should follow the directions for use supplied by the negative pressure wound therapy system manufacturer concerning the protection of exposed underlying structures.

Contraindications

- Avoid using Silverlon® flex Wound Contact Negative Pressure Dressings on patients with known sensitivity to silver or nylon.
- Silverlon® flex Wound Contact Negative Pressure Dressings are not intended for use on 3rd degree burns.

Adverse Reactions

Silverlon® flex Wound Contact Negative Pressure Dressings have been subjected to independent standard in vitro and in vivo biocompatibility tests, including cytotoxicity, sensitization and intracutaneous reactivity. All tests were performed in accordance with the International Standard Organization (ISO) 10993 Standard Series for Biological Evaluation of Medical Devices. The results of the studies indicated that Silverlon® flex Wound Contact Negative Pressure Dressings are safe for their intended use.

Instructions for Use

- Cleanse wound with sterile water, distilled water, or normal saline, removing necrotic debris or eschar as needed per local protocol.
- Select the Silverlon® flex Wound Contact Negative Pressure Dressing size that fully covers the wound up to the wound margins; trim the dressing edges as necessary.
- Activate the Silverlon® flex dressing by thoroughly moistening with sterile water, distilled water, or normal saline.
- Position the Silverlon® flex dressing directly over wound, with either silver side in contact with the skin; secure the dressing in place using a secondary dressing per local protocol.
 - For exudating wounds, use an absorbent secondary dressing of choice.
 - For dry wounds, use a moisture-donating secondary dressing such as hydrocolloid or pre-moistened foam or gauze.
- Periodically check the edges of the Silverlon® flex dressing to ensure that it is maintained in a moist condition.
- For NPWT Use - Position the Silverlon® flex Dressing directly over the wound, in contact with the skin; place foam or gauze fill media of choice on top of the Silverlon® flex dressing and secure in place according to the instructions supplied by the manufacture of the negative pressure wound therapy equipment.
- Silverlon® flex Wound Contact Negative Pressure Dressings may be used for up to 7 days, but may require more frequent changing depending on wound condition and exudate build-up.
 - In conjunction with negative pressure replace the Silverlon® flex dressing whenever the foam / gauze fill media are replaced.
- To remove Silverlon® flex Wound Contact Negative Pressure Dressings, first remove the foam / gauze fill media per local protocol, and then gently depress surrounding skin while lifting the edges of the Silverlon® flex dressing.
 - If sticking of the dressing to the wound occurs, moisten the dressing as needed with sterile water, distilled water, or normal saline, until it can be easily removed by gently lifting the corners.

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