SILVERLON

safe, strong, simple

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SILVERLON®

safe, strong, simple™

Silverlon® Wound Packing Strips (WPS)

EC REP

STERILE EO

Emergo Europe
Molenstraat 15
2513 BH The Hague, The Netherlands
Tel: +(31) (0) 70 345-8570
Fax: +(31) (0) 70 346-7299

Argentum Medical
2571 Kaneville Ct
Geneva, IL 60134 USA
Tel. (+1) 888.551.0188
Fax. (+1) 888.558.9923
www.silverlon.com

LC-CE-IFU-WPS Rev. 03

SILVERLON

safe, strong, simple

Silverlon® Wundstreifen-Verbände**German****SILVERLON**

safe, strong, simple

Silverlon® Wound Packing Strips (WPS)**English**

Device Description

Silverlon® Wound Packing Strips are sterile, single layer, non-adherent, absorbent, antimicrobial dressings, intended for use in wound management, to reduce risk of wound infection.

Available as:

Code	Size
WPS112	2,5 x 14 cm
WPS124	2,5 x 61 cm

Silverlon® Wound Packing Strips contain silver to provide an ancillary, local antimicrobial effect to reduce the risk of wound infection by organisms sensitive to silver. Ionic silver is a broad spectrum antimicrobial, which has been shown to be effective against a range of pathogens associated with wound infection. Silverlon® dressings have been tested in vitro and found effective against microorganisms such as: Enterococcus faecalis, Staphylococcus aureus, Pseudomonas aeruginosa, Acinetobacter baumannii, Enterobacter cloacae, Staphylococcus epidermidis, and Klebsiella pneumoniae.

In vitro testing has shown Silverlon® dressings to be MRI compatible.

Indications

Silverlon® Wound Packing Strips are professional use wound dressings indicated for local management of:

- 1st and 2nd degree burns
- Donor and graft site wounds
- Traumatic and surgical wounds
- Diabetic, pressure and venous ulcers
- Infected wounds*

Silverlon® Wound Packing Strips may be used on those infected wound types listed above. Where used on infected wounds, the infection should be treated as per local clinical protocol.

Silverlon® Wound Packing Strips are indicated for a total duration of contact (i.e. where treatment involves consecutive application of individual dressings) of up to 30 days.

Warnings

- Do not use past expiration date on the product packaging.
- Do not use if pouch is damaged or open.
- Clinicians / Healthcare Professionals should be aware that there are very limited data on prolonged and repeated use of silver containing dressings, particularly in children and neonates.
- Silverlon® dressings should not be used during pregnancy or lactation because of the lack of data concerning the effects of silver on reproduction.
- Do not use petroleum-based ointments or creams under Silverlon® Wound Packing Strips.
- Do not moisten Silverlon® Wound Packing Strips with hydrogen peroxide or povidone iodine.
- Silverlon® Wound Packing Strips are intended for single use only and should not be re-used. Re-use may adversely affect the dressing's performance characteristics and also presents a risk of infection to the patient.
- Total duration of use of Silverlon® Wound Packing Strips (i.e. where treatment involves consecutive application of individual dressings) should not exceed 30 days.

Contraindications

- Avoid using Silverlon® Wound Packing Strips on patients with known sensitivity to silver or nylon.
- Silverlon® Wound Packing Strips are not intended for use on 3rd degree burns.

Adverse Reactions

Silverlon® Wound Packing Strips have been subjected to independent standard in vitro and in vivo biocompatibility tests, including cytotoxicity, sensitization and intracutaneous reactivity. All tests were performed in accordance with the International Standard Organization (ISO) 10993 Standard Series for Biological Evaluation of Medical Devices. The results of these studies indicated that Silverlon® Wound Packing Strips are safe for their intended use.

Instructions for Use

- Cleanse wound with sterile water, distilled water, or normal saline, removing necrotic debris or eschar as needed per local protocol.
- To allow for easy removal select the dressing size that will leave at least ½" of Silverlon® Wound Packing Strip protruding from the wound.
- Trim dressing length as necessary.
- Activate Silverlon® Wound Packing Strip by thoroughly moistening the dressing with sterile water, distilled water, or normal saline.
- Gently pack wound per local protocol ensuring maximum direct contact between Silverlon® Wound Packing Strip and the wound bed.
- Cover and secure Silverlon® Wound Packing Strip with conventional absorbent and non-occlusive dressings per local protocol.
- To remove Silverlon® Wound Packing Strip, first remove the outer dressing per local protocol, then gently grasp and pull the protruding section of Silverlon® Wound Packing Strip.
 - o If sticking of the dressing to the wound occurs, moisten the dressing as needed with sterile water, distilled water, or normal saline, until it can be easily removed by gently pulling the Silverlon® Wound Packing Strip.
- Silverlon® Wound Packing Strips are intended for use for up to 7 days, but dressing change frequency will depend on wound condition and the amount of exudates.
 - o Remove and replace Silverlon® Wound Packing Strip whenever it becomes heavily coated with wound exudates or wound debris.

Compressas Silverlon® Wound Packing Strips (WPS)

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