SILVERLON®

safe, strong, simple

Apósitos tipo tira (WPS) Silverlon®

Español

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SILVERLON®

safe, strong, simple

Silverlon® Wound Packing Strips (WPS)

English

Device Description

Silverlon® Wound Packing Strips are sterile, single layer, non-adherent, absorbent, antimicrobial dressings, intended for use in wound management, to reduce risk of wound infection.

Available as:

Code	Size
WPS112	3 cm x 30 cm
WPS124	3 cm x 61 cm

Silverlon® Wound Packing Strips contain silver to provide an ancillary, local antimicrobial effect to reduce the risk of wound infection by organisms sensitive to silver. Ionic silver is a broad spectrum antimicrobial, which has been shown to be effective against a range of pathogens associated with wound infection. Silverlon® dressings have been tested in vitro and found effective against microorganisms such as: Enterococcus faecalis, Staphylococcus aureus, Pseudomonas aeruginosa, Acinetobacter baumannii, Enterobacter cloacae, Staphylococcus epidermidis, and Klebsiella pneumoniae.

In vitro testing has shown Silverlon® dressings to be MRI compatible.

Indications

Silverlon® Wound Packing Strips are professional use wound dressings indicated for local management of:

- 1st and 2nd degree burns
- Donor and graft site wounds
- Traumatic and surgical wounds
- Diabetic, pressure and venous ulcers
- Infected Wounds*

*Silverlon® Wound Packing Strips may be used on those infected wound types listed above. Where used on infected wounds, the infection should be treated as per local clinical protocol.

Silverlon® Wound Packing Strips are indicated for a total duration of contact (i.e. where treatment involves consecutive application of individual dressings) of up to 30 days.

Warnings

- Do not use past expiration date on the product packaging.
- Do not use if pouch is damaged or open.
- Clinicians / Healthcare Professionals should be aware that there are very limited data on prolonged and repeated use of silver containing dressings, particularly in children and neonates.
- Silverlon® dressings should not be used during pregnancy or lactation because of the lack of data concerning the effects of silver on reproduction.
- Do not use petroleum-based ointments or creams under Silverlon® Wound Packing Strips.
- Do not moisten Silverlon® Wound Packing Strips with hydrogen peroxide or povidone iodine.
- Silverlon® Wound Packing Strips are intended for single use only and should not be re-used. Re-use may adversely affect the dressing's performance characteristics and also presents a risk of infection to the patient.
- Total duration of use of Silverlon® Wound Packing Strips (i.e. where treatment involves consecutive application of individual dressings) should not exceed 30 days.

Contraindications

- Avoid using Silverlon® Wound Packing Strips on patients with known sensitivity to silver or nylon.
- Silverlon® Wound Packing Strips are not intended for use on 3rd degree burns.

Compressas Silverlon® Wound Packing Strips (WPS)

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