SILVERLON[®]
safe, strong, simple™

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English

Device Description

Silverlon® Wound Contact Dressing, Burn Contact Dressing, Elastic Burn Wrap and Digit Sleeve Dressings are sterile, single layer, non-adherent, antimicrobial barrier dressings, intended for topical use in wound management, to reduce risk of wound infection.

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Silverlon® Wound Contact Dressing, Burn Contact Dressing, Elastic Burn Wrap and Digit Silverlon® Wound Contact Dressing and Burn Contact Dressing contain silver to provide an ancillary, local antimicrobial effect to reduce the risk of wound infection by organisms sensitive to silver. Ionic silver is a broad spectrum antimicrobial, which has been shown to be effective against a range of pathogens associated with wound infection. Silverlon® Wound Contact Dressing and Burn Contact Dressing have been tested in vitro and found to be effective against micro organisms such as: Enterococcus faecalis, Staphylococcus aureus, Pseudomonas aeruginosa, Acinetobacter Baumannii, Enterobacter cloacae, Staphylococcus epidermidis, and Klebsiella pneumoniae.

In vitro testing has shown Silverlon® dressings to be MRI-compatible.

Indications

Silverlon® Wound Contact Dressing, Burn Contact Dressing, Elastic Burn Wrap and Digit Sleeve Dressing are professional use wound dressings indicated for local management of:

- 1st and 2nd degree burns
- Donor and graft site wounds
- Traumatic and surgical wounds
- Diabetic, pressure and venous ulcers
- Infected wounds*

*Silverlon® Wound Contact Dressing, Burn Contact Dressing, Elastic Burn Wrap and Digit Sleeve Dressing may be used on those infected wound types listed above. Where used on infected wounds, the infection should be treated as per local clinical protocol.

Silverlon® Wound Contact Dressing, Burn Contact Dressing, Elastic Burn Wrap and Digit Sleeve Dressing are indicated for a total duration of contact (i.e. where treatment involves consecutive application of individual dressings) of up to 30 days.

Warnings

- Do not use past expiration date on the product packaging.
- Do not use if pouch is damaged or open.
- Clinicians / Healthcare Professionals should be aware that there are very limited data on prolonged and repeated use of silver-containing dressings, particularly in children and neonates.
- Silverlon® Wound Contact Dressing, Burn Contact Dressing, Elastic Burn Wrap and Digit Sleeve Dressing should not be used during pregnancy or lactation because of the lack of data concerning the effects of silver on reproduction.
- Do not use petroleum-based ointments or creams under Silverlon® Wound Contact Dressing, Burn Contact Dressing, Elastic Burn Wrap and Digit Sleeve Dressing
- Do not moisten Silverlon® Wound Contact Dressing, Burn Contact Dressing, Elastic Burn Wrap and Digit Sleeve Dressing with hydrogen peroxide or povidone iodine.
- Zum Entfernen des Silverlon® Wound Contact Dressing, Elastic Burn Wrap and Digit Sleeve Dressing are intended for single use only and should not be reused. Re-use may adversely affect the dressing’s performance characteristics and also presents a risk of infection to the patient.

• Total duration of use of Silverlon® Wound Contact Dressing, Burn Contact Dressing, Elastic Burn Wrap and Digit Sleeve Dressing (i.e. where treatment involves consecutive application of individual dressings) should not exceed 30 days.

Contraindications

• Avoid using Silverlon® Wound Contact Dressing, Burn Contact Dressing, Elastic Burn Wrap and Digit Sleeve Dressing on patients with known sensitivity to silver or nylon.

• Silverlon® Wound Contact Dressing, Burn Contact Dressing, Elastic Burn Wrap and Digit Sleeve Dressing are not intended for use on 3rd degree burns.

Adverse Reactions

Silverlon® Wound Contact Dressing, Burn Contact Dressing, Elastic Burn Wrap and Digit Sleeve Dressing have been subjected to independent standard in vitro and in vivo biocompatibility tests, including cytotoxicity, sensitization, and intracutaneous reactivity. All tests were performed in accordance with the International Standard Organization (ISO) 10993 Standard Series for Biological Evaluation of Medical Devices. The results of the studies indicated that Silverlon® Wound Contact Dressing, Burn Contact Dressing, Elastic Burn Wrap and Digit Sleeve Dressing are safe for their intended use.

Instructions for Use

- Cleanse wound with sterile water, distilled water, or normal saline removing necrotic debris or eschar as needed per local protocol.
- Select the dressing size that overlaps the wound margins by 1–2 cm.
- Activate the Silverlon® dressing by thoroughly moistening with sterile water, distilled water or normal saline.
- Position the Silverlon® dressing directly over wound, with either silver side in contact with the skin; secure the dressing in place using a secondary dressing per local protocol.
 - For exudating wounds, use an absorbent secondary dressing of choice.
 - For dry wounds, use a moisture donating secondary dressing such as hydrocolloid or pre-moistened foam or gauze.
- Periodically check the edges of the Silverlon® dressing to ensure that it is maintained in a moist condition.
- Silverlon® Wound Contact Dressing, Burn Contact Dressing, Elastic Burn Wrap and Digit Sleeve Dressing may be used for up to 7 days, but may require more frequent changing depending on wound condition and exudate buildup.
- To remove Silverlon® Wound Contact Dressing, Burn Contact Dressing, Elastic Burn Wrap and Digit Sleeve Dressing, first remove the outer secondary dressing per local protocol, and then gently depress surrounding skin while lifting the dressing edges.
 - If sticking of the dressing to the wound occurs, moisten the dressing as needed with sterile water, distilled water or normal saline until it can be easily removed by gently lifting the corners.

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