**SILVERLON**  
safe, strong, simple™

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**English**

#### Device Description

Silverlon® Wound Pad and Burn Pad Dressings are sterile, multi-layer, non-adherent, absorbent, antimicrobial barrier dressings, intended for topical use in wound management, to reduce risk of wound infection.

**Available as:**

Code	Pad Size	Code	Pad Size	Code	Pad Size
WPD22	5 x 5 cm	WPD212	5 x 30 cm	BPD44	10 x 10 cm
WPD23	5 x 7 cm	WPD38	7 x 20 cm	BPD48	10 x 20 cm
WPD26	5 x 15 cm	WPD310	7 x 25 cm	BPD816	20 x 40 cm
WPD28	5 x 20 cm	WPD316	7 x 40 cm		
WPD210	5 x 25 cm	WPD44	10 x 10 cm		

Silverlon® Acute Burn Glove Dressing contains silver to provide an ancillary, local antimicrobial effect to reduce the risk of wound infection by organisms sensitive to silver. Ionic silver is a broad spectrum antimicrobial, which has been shown to be effective against a range of pathogens associated with wound infection. Silverlon® dressings have been tested in vitro and found effective against micro organisms such as: Enterococcus faecalis, Staphylococcus aureus, Pseudomonas aeruginosa, Acinetobacter Baumannii, Enterobacter cloacae, Staphylococcus epidermidis and Klebsiella pneumoniae.

In vitro testing has shown Silverlon® dressings to be MRI compatible.

#### Indications

Silverlon® Wound Pad Dressing and Burn Pad Dressings are professional use wound dressings intended for management of:

- Trauma wounds
- Surgical wounds
- Donor and graft site wounds
- 1st and 2nd degree burns
- Chronic wounds (diabetic, pressure and venous ulcers)

Silverlon® Wound Pad and Burn Pad Dressings may be used on infected wounds. Where used on infected wounds, the infection should be treated as per local clinical protocol. Silverlon® Wound Pad and Burn Pad Dressings Dressings are indicated for a total duration of contact (i.e. where treatment involves consecutive application of individual dressings) of up to 30 days

#### Warnings

- Do not use past expiration date on the product packaging.
- Do not use if pouch is damaged or open.
- Clinicians / Healthcare Professionals should be aware that there are very limited data on prolonged and repeated use of silver-containing dressings, particularly in children and neonates.
- Silverlon® dressings should not be used during pregnancy or lactation because of the lack of data concerning the effects of silver on reproduction.
- Do not use petroleum-based ointments or creams under Silverlon® Wound Pad and Burn Pad Dressing.
- Do not moisten Silverlon® Wound Pad and Burn Pad Dressings with hydrogen peroxide, or povidone iodine.
- Silverlon® Wound Pad and Burn Pad Dressings are intended for single use only and should not be re-used. Re-use may adversely affect the dressing's performance characteristics and also presents a risk of infection to the patient.
- Total duration of use of Silverlon® Wound Pad and Burn Pad Dressings (i.e. where treatment involves consecutive application of individual dressings) should not exceed 30 days.

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