Opatrunki Silverlon® Wound Pad Dressings Opatrunki Silverlon® Burn Pad Dressings

Polskie

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English

Device Description

Silverlon® Wound Pad and Burn Pad Dressings are sterile, multi-layer, non-adherent, absorbent, antimicrobial barrier dressings, intended for topical use in wound management, to reduce risk of wound infection.

Available as:

Code	Pad Size	Code	Pad Size	Code	Pad Size
WPD22	5 x 5 cm	WPD212	5 x 30 cm	BPD44	10 x 12 cm
WPD23	5 x 8 cm	WPD38	8 x 20 cm	BPD48	10 x 20 cm
WPD26	5 x 15 cm	WPD310	8 x 25 cm	BPD816	20 x 40 cm
WPD28	5 x 20 cm	WPD316	8 x 40 cm		
WPD210	5 x 25 cm	WPD44	10 x 12 cm		

Silverlon® Wound Pad and Burn Pad Dressing contain silver to provide an ancillary, local antimicrobial effect to reduce the risk of wound infection by organisms sensitive to silver. Ionic silver is a broad spectrum antimicrobial, which has been shown to be effective against a range of pathogens associated with wound infection. Silverlon® Wound Pad and Burn Pad Dressings have been tested in vitro and found effective against microorganisms such as: *Enterococcus faecalis*, *Staphylococcus aureus*, *Pseudomonas aeruginosa*, *Acinetobacter Baumannii*, *Enterobacter cloacae*, *Staphylococcus epidermidis* and *Klebsiella pneumoniae*.

In vitro testing has shown Silverlon® dressings to be MRI compatible.

Indications

Silverlon® Wound Pad Dressing and Burn Pad Dressings are professional use wound dressings intended for management of:

- 1st and 2nd degree burns
- Donor and graft site wounds
- Diabetic, pressure and venous ulcers
- Trauma and surgical wounds
- Infected wounds*

*Silverlon® Wound Pad and Burn Pad Dressings may be used on the infected wounds listed above. Where used on infected wounds, the infection should be treated as per local clinical protocol.

Silverlon® Wound Pad and Burn Pad Dressings Dressings are indicated for a total duration of contact (i.e. where treatment involves consecutive application of individual dressings) of up to 30 days.

Warnings

- Do not use past expiration date on the product packaging.
- Do not use if pouch is damaged or open.
- Clinicians / Healthcare Professionals should be aware that there are very limited data on prolonged and repeated use of silver-containing dressings, particularly in children and neonates.
- Silverlon® dressings should not be used during pregnancy or lactation because of the lack of data concerning the effects of silver on reproduction.
- Do not use petroleum-based ointments or creams under Silverlon® Wound Pad and Burn Pad Dressing.
- Do not moisten Silverlon® Wound Pad and Burn Pad Dressings with hydrogen peroxide, or povidone iodine.
- Silverlon® Wound Pad and Burn Pad Dressings are intended for single use only and should not be re-used. Re-use may adversely affect the dressing's performance characteristics and also presents a risk of