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### Silverlon® Wound Packing Strips (WPS)

### English

<b>SILVERLON</b>
<i>PERFORMANCE. POWER. SPEED.</i>

**Device Description**

Silverlon® Wound Packing Strips are sterile, single layer, non-adherent, absorbent, antimicrobial dressings, intended for use in wound management, to reduce risk of wound infection.

**Available as:**

<b>Code</b>	<b>Size</b>
WPS112	3 cm x 30 cm
WPS124	3 cm x 61 cm

Silverlon® Wound Packing Strips contain silver to provide an ancillary, local antimicrobial effect to reduce the risk of wound infection by organisms sensitive to silver. Ionic silver is a broad spectrum antimicrobial, which has been shown to be effective against a range of pathogens associated with wound infection. Silverlon® dressings have been tested in vitro and found effective against microorganisms such as: Enterococcus faecalis, Staphylococcus aureus, Pseudomonas aeruginosa, Acinetobacter baumannii, Enterobacter cloacae, Staphylococcus epidermidis, and Klebsiella pneumoniae.

In vitro testing has shown Silverlon® dressings to be MRI compatible.

**Indications**

Silverlon® Wound Packing Strips are professional use wound dressings indicated for local management of:

- 1st and 2nd degree burns
- Donor and graft site wounds
- Traumatic and surgical wounds
- Diabetic, pressure and venous ulcers
- Infected Wounds\*

\*Silverlon® Wound Packing Strips may be used on those infected wound types listed above. Where used on infected wounds, the infection should be treated as per local clinical protocol.

Silverlon® Wound Packing Strips are indicated for a total duration of contact (i.e. where treatment involves consecutive application of individual dressings) of up to 30 days.

**Warnings**

- Do not use past expiration date on the product packaging.
- Do not use if pouch is damaged or open.
- Clinicians / Healthcare Professionals should be aware that there are very limited data on prolonged and repeated use of silver containing dressings, particularly in children and neonates.
- Silverlon® dressings should not be used during pregnancy or lactation because of the lack of data concerning the effects of silver on reproduction.
- Do not use petroleum-based ointments or creams under Silverlon® Wound Packing Strips.
- Do not moisten Silverlon® Wound Packing Strips with hydrogen peroxide or povidone iodine.
- Silverlon® Wound Packing Strips are intended for

*PERFORMANCE. POWER. SPEED.*

# Silverlon® Wound Packing Strips (WPS)

					
					
<b>Emergo Europe, Westervoortsedijk 60; 6827 AT Arnhem The Netherlands Tel. +(31) (0) 70 345-8570</b>					<b>Argentum Medical, LLC 2571 Kaneville Ct Geneva, IL 60134 USA Tel. (+1) 888.551.0188 Fax. (+1) 888.558.9923 www.silverlon.com</b>
<b>LC-CE-IFU-WPS-03; Rev B</b>					

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